FLEX
Monoclonal Mouse
Anti-Human
Mast Cell Tryptase
Clone AA1
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

Code IS640

ENGLISH	
Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human Mast Cell Tryptase, Clone AA1, Ready-to-Use (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody labels human mast cells and is useful for the identification of very atypical or immature mast cells (MC) in mast cell leukemia, and for the detection of small, even minute, dense focal MC infiltrates in staging procedures in patients with known cutaneous mastocytosis (1). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Human mast cell tryptases (EC 3.4.21.59) comprise a family of trypsin-like neutral serine proteases that are predominantly expressed in mast cells (2). In its enzymatically active form, mast cell tryptase exists as a non-covalently linked tetramer of 132 kDa (2, 3). Mast cell tryptase is capable of degrading vasoactive intestinal peptide and activating prekallikrein as well as generating kinins, all important mediators involved in bronchoconstriction and airway hyperresponsiveness, which are major contributors to allergic airway disease (2). Mast cells are activated by a number of stimuli, including antigen, superoxides, complement protein, neuropeptides and lipoproteins, resulting in activation and degranulation. Mast cells play an active role in such diverse diseases as atherosclerosis, asthma, arthritis, bile duct fibrosis, malignancy and pulmonary fibrosis (4, 5). Identification of mast cells through staining of tissues with antibodies specific for human mast cell tryptase has been useful in identification of focal and diffuse MC infiltrates in primary MC disease and mastocytosis. Demonstration of mast cell tryptase is essential for the identification of highly atypical, hypogranulated or even non-metachromatic MCs, especially in MC leukemia, and for the detection of small or minute MC infiltrates (1). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide. <u>Clone</u> : AA1 (3). <u>Isotype</u> : IgG1, kappa.
Immunogen	Mast cell tryptase isolated from human lung (3).
Specificity	In Western blotting of a tryptase-rich extract of human lung under reducing conditions, the antibody labels a band of approximately 32.5 kDa corresponding to human mast cell tryptase (3). In indirect ELISA the antibody reacts with purified human mast cell tryptase (3). Antibody-binding to the mast cell tryptase does not affect enzyme activity (3). Periodate treatment of mast cell tryptase has no effect on antibody labeling, indicating that carbohydrate epitopes are not involved in antibody recognition. Nor is antibody recognition affected by the presence of heparin (3).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections</u> : Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections</u> : Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium.
Staining procedure including materials required but not supplied	The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: Trypt (without counterstaining) or TryptH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument.

(117438-002)

IS640/EFG/MNI/2009.12.04 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

Staining interpretation
Performance characteristics

If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services.
Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics".
Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin, (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended.
Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include tonsil and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse, (Dako Autostainer/Autostainer Plus) (Code IS750).

Cells labeled by the antibody display granular cytoplasmic staining (6).

Normal tissues: The antibody specifically labels mast cells (MC) in a number of human tissues, including lung, tonsil, colon, gastric mucosa and skin, whereas only a few MCs are detectable in the pituitary gland. The antibody does not label basophils or any other normal blood leukocytes, likewise, no labeling has been observed in keratinocytes (6). Mast cells in the tonsil show a moderate to strong staining reaction.
Abnormal tissues: In a study of more than 150 cases of different MC disorders, it was demonstrated that <4 MCs/mm² were labeled in bone marrow in most normal and reactive states, whereas >5 but <100 were labeled in myelodysplastic syndromes. More than 100 MCs/mm² were detected in MC neoplasias, and cases with a much higher number (maximum 2655 MCs/mm²) were often observed (1).

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- Horny H-P, Valent P. Diagnosis of mastocytosis: general histopathological aspects, morphological criteria, and immunohistochemical findings. Leuk Res 2001;25:543-51.
- Abraham WM. Tryptase: potential role in airway inflammation and remodelling. Am J Physiol Lung Cell Mol Physiol 2002;282:L193-L196.
- Walls AF, Bennett AR, McBride HM, Glennie MJ, Holgate ST, Church MK. Production and characterization of monoclonal antibodies specific for human mast cell tryptase. Clin Exp Allergy 1990;20:581-9.
- Krishnaswamy G, Kelley J, Johnson D, Youngberg D, Stone W, Huang S-K, et al. The human mast cell: functions in physiology and disease. Front Biosci 2001;6:D1109-27.
- Koda W, Harada K, Tsuneyama K, Kono N, Sasaki M, Matsui O, et al. Evidence of the participation of peribiliary mast cells in regulation of the peribiliary vascular plexus along the intrahepatic biliary tree. Lab Invest 2000;80:1007-17.
- Walls AF, Jones DB, Williams JH, Church MK, Holgate ST. Immunohistochemical identification of mast cells in formaldehyde-fixed tissue using monoclonal antibodies specific for tryptase. J Pathol 1990;162:119-26.

Explanation of symbols/ Explication des symboles/ Erläuterung der Symbole			
REF	Catalogue number Référence catalogue Bestellnummer	2 °C  8 °C Temperature limitation Limites de température Zulässiger Temperaturbereich	 Use by Utiliser avant Verwendbar bis
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Contains sufficient for <n> tests Contenu suffisant pour <n> tests Inhalt ausreichend für <n> Tests	 Manufacturer Fabricant Hersteller
	Consult instructions for use Voir les instructions d'utilisation Gebrauchsanweisung beachten	LOT Batch code Numéro de lot Chargenbezeichnung	

(117438-002)

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

IS640/EFG/MNI/2009.12.04 p. 4/4